




Germany

# RAMAN & WEIL PVT LTD

*Since 1989*

**Advanced Disinfectants for  
Instruments & Endoscopes  
Product Range**

**J.S. ANAND  
GENERAL MANAGER  
RAMAN & WEIL PVT. LTD.**





# Excellence In Infection Control



31 years of association with **Bode Chemie, Germany**, R&W has more than **44 products**, covering the entire spectrum of Hi Tech Disinfection and Antiseptic needs of all Medical spheres: Hospitals, Nursing Homes, Dental Clinics, Pathological laboratories & Pharmaceutical Sterile Sections.





# **Excellence in Infection Control**

**In technical collaboration with**  
**BODE CHEMIE, Germany**



- **DIN ISO 13485**
- **DIN ISO 9001**
- **DIN ISO 14001**
- **GMP Certificate**
- **EC Certificate of Conformity**
- **Manufacturer's Authorization**
- **Certificate Environmental Partnership Program**
- **WHO GMP Certificate in all 3 units**





The number 1 trusted name in  
Hand Disinfection :

**Sterillium**

The International Gold Std.







RAMAN & WEIL Pvt. Ltd.

EXCELLENCE IN INFECTION CONTROL



# Sterillium

## The International Gold Std.

Passes EN 1500, EN 12791,  
EN 13727, EN 13624, EN 14348,  
EN 1040, EN 1275  
Brand of the Century, Europe

**International Gold Standard for  
alcohol hand rubs**

Approved & used in 57 countries

**As per World Health Organization,  
recommendation**

2016 – 4 Billion hand disinfection  
procedures conducted, 9000 tonnes

**0.2 micron filtration conducted  
under GMP conditions**

Manufactured as per European  
Pharmacopeia, highest purity of  
alcohol

**Skin Hydration increases up to 30%**


Highly efficient formula







# Range of Effective Instrument & Endoscope Disinfectants





# KORSOLEX PLUS

- Aldehyde free instrument disinfectant (Conc.)
- Ideal for heat sensitive and heat resistance instruments
- With powerful cleaning effect, corrosion inhibitors, pH regulators and fragrance.
- Very safe, as this product, will not disturb the lens cement which holds the lenses, so the focus remains accurate, in spite of frequent use.
- STANDARD SOLUTION- Use 2% only
- STABILITY  
Diluted solution can be used up to 7 days, including high protein loads  
Shelf life of the concentrate is 3 yrs.

Excellence in  
Infection Control



## Korsolex® Plus

Aldehyde - free instrument disinfectant for heat-sensitive and heat-resistant instruments. Powerful cleaning effect.



Helicobacter Pylori  
(Rem magnification)



# KORSOLEX PLUS

## CHEMICAL COMPOSITION

- 5TH Generation QAC
- With very powerful cleanser
- Manufactured as per European Pharmacopeia
- Listed in the German Association for Applied Hygiene
- CE Labeling in accordance with the Medical Device Directive.

## MATERIAL COMPATIBILITY

- Excellent material compatibility on glass, ceramic, stainless steel, non-ferrous metals, aluminum, hard rubber, hard plastic, silicone, Makrolone, & Plexiglas.

### Korsolex® Plus

Aldehyde - free instrument disinfectant for heat-sensitive and heat-resistant instruments. Powerful cleaning effect.

#### Proven efficacy

Bacteria and Fungi		5 min	15 min
EN Phase 2 / Step 2 Efficacy according to EN Phase 2 / Step 2 (practical tests) tested under clean / dirty conditions	Bactericidal (EN 14561) - dirty conditions Yeasticidal (EN 14562) - dirty conditions	0.1%	0.25%
EN Phase 2 / Step 1 Efficacy according to EN Phase 2 / Step 1 (suspension tests) tested under clean / dirty conditions	Bactericidal (EN 13727) - dirty conditions Yeasticidal (EN 13624) - dirty conditions	0.25%	0.1%
EN Phase 1 / Basic test Appraised efficacy according to EN Phase 1 (Basic tests / suspension tests) without contamination does not define the applicability of a product for a specific purpose	Bactericidal (EN 1040) Yeasticidal (EN 1275) Fungicidal (EN 1275)	0.05% 0.1%	0.5%
VAH Certified application recommendation for prophylactic instrument disinfection from the Association for Applied Hygiene (VAH); Based on suspension and practical tests; tested under clean (i.e. cleaned instruments) dirty conditions (i.e. used / dirty instruments)	Bactericidal / Yeasticidal / Tuberculocidal - dirty conditions		30 mins 2.0%
DGHM Appraised efficacy against bacteria (in accordance with the German Society for Hygiene and Microbiology [DGHM]); within the certified bactericidal efficacy	+ Antibiotic resistant bacteria + Helicobacteria pylori	1 min 1%	15 mins 0.5%
<b>Viruses</b>			
Efficacy against viruses (acc. to German Society for the Control of Viral Diseases - DVV)	Virucidal against enveloped viruses (incl. HBV, HIV, HCV)	5 mins	2%
Appraised efficacy against enveloped viruses (DVV)	BVDV Vaccinia virus	2% 2%	
Appraised efficacy against enveloped viruses (in accordance with DVV)	Herpes - simplex virus (Type 1)	2%	
Appraised efficacy against non-enveloped viruses (in accordance with DVV)	Rota virus	2%	

#### Composition

100 g concentrate contain:  
Dodecylbispropylene triamine - 9.2 g  
Didecylidimethyl ammonium Chloride 13.0 g  
Other ingredients : Surfactants, corrosion inhibitors, foam regulators, pH regulators, fragrance.

#### Microbiology

Bactericidal / Yeasticidal / Tuberculocidal / Virucidal against enveloped viruses (incl. HBV, HIV, HCV) & non-enveloped viruses (Rota)

#### Areas of application

Korsolex® plus can be used for the cleaning disinfection of heat-resistant and heat-sensitive instruments in the manual treatment, circulation procedures as well as in the ultrasonic bath. Safe for cleaning disinfection of CO<sub>2</sub> incubators, Laminar flow hoods, Ultrasound probes, etc. present in IVF Clinics and Laboratories.

#### Directions for use

**Immersion bath procedure**  
Korsolex® plus is supplied as concentrate. Open or take instruments and devices apart as far as possible. Make sure the instruments are completely covered by Korsorex plus solution without any hollows or bubbles.

#### Product properties

Virucidal against enveloped viruses (incl. HBV, HIV, HCV) exceptional cleaning power high material compatibility pleasant smell suitable for ultrasonic cleaning and disinfection

#### Listing

Certification / list issued by the German Association for Applied Hygiene (VAH) CE labelling in accordance with the Medical Device Directive (MDD)

#### Stability

- **After opening**  
in tightly closed container : 12 months
- **Stability of fresh working solution**  
in closed containers : 4 weeks
- **Standing time of the used working solution :**  
The use-solution can be utilized for 7 days of including low to high protein loads. This is also applicable for ultrasonic baths.

#### Endorsed by



#### Ultrasonic bath

Korsolex® plus can be used in all stainless steel ultrasonic baths. Follow the sonication time stated in the instruments manufacturer's instructions and make sure the solution is not heated, the exposure time has to be followed.

#### Rinse phase

After disinfection, adequate rinsing shall be carried out. In practice, rinse times of 10 minutes have proved most appropriate. The water should be of at least drinking water quality.

#### Compatibility

The material compatibility of the instrument disinfectant extends to a wide range of materials, e.g. glass, ceramic, stainless steel, non-ferrous metals, aluminium, hard rubber and hard plastic, rubber (including latex and silicone) and also plastics (including Makrolon and Plexiglas).

**Standard Solution : Use 2% only**  
Always prepare the solution with max. room temperature water

Prepared Solution	Korsolex® Plus	Potable Drinking Water
5L	100 ml	4.900 L
2L	40 ml	1.960 L
1L	20 ml	0.980 L

#### Packing :

100 ml : 50 bottles / Carton



#### Organic loads for testing instrument disinfectants

##### Test conditions - instrument disinfection :

**Clean conditions (following EN standards) :**  
Products tested under clean conditions only are recommended for pre-cleaned instruments and/or instruments that are minimally contaminated with organic and / or inorganic substances.

**Dirty conditions (following EN standards) :**  
In case you cannot rule out contamination with blood, only use products with proven activity under dirty conditions.

Exclusively Manufactured & Marketed by:

**RAMAN & WEIL PVT. LTD.**

Registered Office : 15, Chateau Marine, B-Rd., Churchgate, Mumbai - 20.  
Administrative Office : Room No. 36/37, 3rd Floor, Mittal Chamber, Opp. Inox Theatre, Nariman Point, Mumbai - 400 021, India  
Phone : 22049527, 22856397 Email : sales@ramanweil.com  
Visit : www.ramanweil.com



In technical collaboration with:

**BODE Chemie GmbH**

Melanchthonstraße 27  
D-22525 Hamburg, Germany  
Must Visit : www.bode-science-center.com  
Recipient : ISO 9001, ISO 14001, ISO 13485  
EC Certificate of Conformity, WHO - GMP

Endorsed by

**FUJINON**  
FUJIFILM

**OLYMPUS**

**STORZ**  
KARL STORZ — ENDOSKOPE

**RICHARD WOLF**





# KORSOLEX PLUS



- EN 1040 – Bactericidal, including antibiotic resistant bacteria
- EN 1275 – Yeasticidal & Fungicidal
- EN 13727 – Bactericidal (incl. diff. MRSA strains)
- EN 13624 – Yeasticidal & Fungicidal
- EN 14561 – Bactericidal, including **Helicobacter Pylori**
- EN 14562 – Yeasticidal
- EN 14348 – Mycobactericidal & Tuberculocidal
- EN 14563 – Mycobactericidal & Tuberculocidal
- EN 14476 – Viricidal (Enveloped & Non Enveloped including Coronavirus)



100ml

*ALL EFFICACY REPORTS, MATERIAL COMPATIBILITY REPORTS, ENDORSEMENTS AVAILABLE ON REQUEST.*



# KORSOLEX RAPID

- Aldehyde based sterilizer concentrate
- For Thermolabile and Thermostable instruments
- Rust Inhibitors and High-Tech Cleansers
- Very safe, as this product, will not disturb the lens cement which holds the lenses, so the focus remains accurate, in spite of frequent use.
- **STANDARD SOLUTION- Use 5% only**

	STRENGTH OF SOLUTION	IMMERSION TIME
STD./ EMERGENCY DISINFECTION (under clean conditions)	5%	5 mins only
STD./ EMERGENCY DISINFECTION (under dirty conditions)	5%	10-15 mins only
TOTAL STERILIZATION	5%	4 hrs


- **STABILITY**  
Diluted solution can be used up to 7 days, including high protein loads.  
Shelf life of the concentrate is 3 yrs.

Excellence in Infection Control



## Korsolex<sup>®</sup> Rapid

Aldehyde-based sterilizer, for heat-sensitive and heat-resistant instruments. Virucidal & Sporicidal activity. Economical concentration.





# KORSOLEX RAPID

## CHEMICAL COMPOSITION

- Complex combination of 2 types of Aldehydes based formulation Concentrate
- With very powerful cleanser
- Manufactured as per European Pharmacopeia
- Listed in the German Association for Applied Hygiene
- CE Labeling in accordance with the Medical Device Directive.

## MATERIAL COMPATIBILITY

- Excellent material compatibility on glass, ceramic, stainless steel, non-ferrous metals, aluminum, hard rubber, hard plastic, silicone, Makrolone, & Plexiglas

Endorsed by



PENTAX



## Korsolex® Rapid

Aldehyde-based disinfectant for heat-sensitive and heat-resistant instruments.

### Proven efficacy

Bacteria and Fungi		
<b>EN Phase 2 / Step 2</b> Efficacy according to EN Phase 2 / Step 2 (Practical tests) tested under clean / dirty conditions	Bactericidal (EN 14561) - clean conditions Yeasticidal (EN 14562) - clean conditions Tuberculocidal - clean conditions Mycobactericidal (EN 14563) - clean conditions	4% - 5 mins     4% - 15 mins
<b>EN Phase 2 / Step 1</b> Efficacy according to EN Phase 2 / Step 1 (Suspension tests) tested under clean / dirty conditions	Bactericidal (EN 13727) - clean conditions Yeasticidal (EN 13624) - clean conditions Fungicidal (EN 13624) - clean conditions Tuberculocidal (EN 14348) - clean conditions Mycobactericidal (EN 14348) - clean conditions	3% - 5 mins     5% - 15 mins
<b>EN Phase 1 / Basic test</b> Appraised efficacy according to EN Phase 1 (basic tests / suspension tests) without contamination does not define the applicability of a product for a specific purpose	Bactericidal (EN 1040) Yeasticidal (EN 1275) Fungicidal (EN 1275)	1% - 1 mins 2% - 1 mins 2% - 5 mins
<b>VAH</b> Certified application recommendation for prophylactic instrument disinfection from the Association for Applied Hygiene (VAH). Based on suspension and practical tests; tested under clean (i.e. cleaned instruments) / dirty conditions (i.e. used / dirty instruments)	Bactericidal / Yeasticidal - clean conditions Fungicidal - clean conditions Tuberculocidal - clean conditions Mycobactericidal - clean conditions	5% - 15 mins
<b>RKI</b> Recognised substance for decontamination acc. to §18 IfSG (Robert Koch-Institut-RKI) Area B - see "Viruses"	Area A - vegetative bacteria; incl. Mycobacteria, fungi and fungi-spores	4% - 15 mins
<b>Bacterial spores</b> Efficacy against bacterial spores	Sporicidal B. subtilis	5% - 4 hrs
<b>Viruses</b> EN Phase 2 / Step 1 Appraised efficacy according to EN Phase 2/Step 1 (suspension tests) tested under clean / dirty conditions	Virucidal (EN 14476) - clean conditions Adenovirus (EN 14476) - clean conditions	1% - 2% 5 mins-15mins
Efficacy against viruses (acc. to German Society for the Control of Viral Diseases -DVV)	Virucidal against enveloped viruses (incl. HBV, HIV, HCV) Virucidal	
Appraised efficacy against non-enveloped viruses (DVV) Appraised efficacy against non-enveloped viruses (in accordance with DVV)	Adenovirus Polyomavirus  Rotavirus	
Recognised substance for decontamination acc. to §18 IfSG (RKI) Area A-see "Bacteria and Fungi"	Area B - enveloped and non-enveloped viruses	3% - 15 min

### Composition

Each 100 gms contain :  
Glutaraldehyde : 15.2 g  
1, 6-Dihydroxy 2, 5-Dioxahexane:19.7g  
Rust Inhibitors and Hi-Tech Cleansors  
4% solution : pH-value approx. 7

### Microbiology

bactericidal (including *Helicobacteria Pylori* EHEC / MRSA), yeasticidal, fungicidal tuberculocidal, mycobactericidal sporicidal, virucidal against enveloped viruses (incl. HBV, HIV, HCV), virucidal

### Areas of Application

**Korsolex Rapid®** displays excellent material compatibility with heat-resistant and heat-sensitive instruments, and is therefore a product of choice in endoscopy. **Korsolex Rapid®** is suitable for all common manual and semi-automatic (cold disinfection) procedure. **Korsolex Rapid®** is supplied as a concentrate. also suitable for **Immersion bath procedure, Circulation process, Ultrasonic bath Chemo-thermal reprocessing of endoscopes**

### Shelf Life:

3 years in conc. form :  
prepared solution 7 days

### Dosage

	Strength of Solution	Immersion Time
Std. / Emergency Disinfection (under clean conditions)	5%	5mins only
Std. / Emergency Disinfection (under dirty conditions)	5%	10-15mins only
Total Sterilization	5%	4 hrs

### Packing

500 ml : 20 bottles / Carton



Endorsed by **STORZ** KARL STORZ - ENDOSKOPE **PENTAX** **RICHARD WOLF**

Supported by Comprehensive proofs of efficacy and scientific - based research and development, Our hygiene and disinfection products ensure best quality.  
[www.bode-science-center.com](http://www.bode-science-center.com)



Exclusively Manufactured & Marketed by:  
**RAMAN & WEIL PVT. LTD.**  
Registered Office : 15, Chateau Marine, B-Rd., Churchgate, Mumbai - 20.  
Administrative Office : Room No. 36/37, 3rd Floor, Mittal Chamber, Opp. Inox Theatre, Nariman Point, Mumbai - 400 021. India  
Phone : 22049527, 22856397 Email : sales@ramanweil.com  
Visit : www.ramanweil.com



In technical collaboration with:  
**BODE Chemie GmbH**  
Melanchthonstraße 27  
D-22525 Hamburg, Germany  
Must Visit : [www.bode-science-center.com](http://www.bode-science-center.com)  
[www.sterillium.com](http://www.sterillium.com)  
Recipient : ISO 9001, ISO 14001, ISO 13485  
EC Certificate of Conformity, WHO - GMP





# KORSOLEX RAPID



- EN 1040 – Bactericidal, including antibiotic resistant bacteria
- EN 1275 – Yeasticidal & Fungicidal
- EN 13727 – Bactericidal (incl. diff. MRSA strains)
- EN 13624 – Yeasticidal & Fungicidal
- EN 14561 – Bactericidal, including **Helicobacter Pylori**
- EN 14562 – Yeasticidal
- EN 14347 – Sporicidal
- EN 14348 – Mycobactericidal & Tuberculocidal
- EN 14563 – Mycobactericidal & Tuberculocidal
- EN 14476 – Viricidal (Enveloped & Non Enveloped  
including Coronavirus)



500ml

*ALL EFFICACY REPORTS, MATERIAL COMPATIBILITY REPORTS, ENDORSEMENTS AVAILABLE ON REQUEST.*



# BODEDEX FORTE

- pH-neutral cleaner for heat sensitive and heat resistant instruments.
- Patented self-acting liquid cleaner
- Reliable and gentle in removal of stubborn residues of blood, proteins, secretions, grease, biofilms and X-ray contrast media.
- Very safe, as this product, will not disturb the lens cement which holds the lenses, so the focus remains accurate, in spite of frequent use.

## PRODUCT PROPERTIES

- Contains special complexing agents to remove water hardness
- Ideal for pre-cleaning and cleaning of instruments having narrow lumen (specialized endoscopes)/ **Duodenoscopes that are heavily contaminated with residues of contrast media.**
- Can be used in Semi-automatic circulation machine/ Ultrasonic Bath/ Removal of contrast media residues
- Used for cleaning of anesthetic instruments like masks, tubes, etc.

Excellence in  
Infection Control



## Bodedex<sup>®</sup> forte

pH-neutral cleaner for heat sensitive and heat resistant instruments.  
Patented self-acting liquid cleaner



### Disinfecting cleaning of instruments



Contaminated instruments and endoscopes put healthcare personnel at risk and should be introduced to the cleaning process as quick as possible. Using a cleaning disinfectant – active against bacteria, fungi and enveloped viruses, e.g. HBV, HIV and HCV – can reduce the risk of infection during wet disposal or manual pre-cleaning.

Combinations with aldehydes or products based on peracetic acid should not be used for pre-cleaning as they fix proteins (1, 2, 3).

An alternative are products based on quaternary ammonium compounds, which are – provided that no aldehyde is added – compatible with all other active agents used in instrument disinfection today, independent of whether the instruments are reprocessed manually or automatically afterwards (4).

1 Anforderungen an die Hygiene bei der Aufbereitung von Medizinprodukten. Empfehlung der Kommission für Krankenhaushygiene und Infektionsprävention beim Robert Koch-Institut (RKI) und des Bundesinstitutes für Arzneimittel und Medizinprodukte (BfArM). Bundesgesundheitsbl. - Gesundheitsforsch. - Gesundheitsschutz, 2001, 44:1115–1126.

2 Arbeitskreis Instrumenten-Aufbereitung. Instrumenten-Aufbereitung richtig gemacht. 9. Ausgabe, 2009, 22-25.

3 Fünf Jahre Empfehlungen der Kommission für Krankenhaushygiene zur Aufbereitung flexibler Endoskope. Bundesgesundheitsbl. 51 (2008):211-220.

4 Zum Beispiel das aldehydfreie Instrumenten-Desinfektionsmittel Bomix plus auf Basis quartärer Ammoniumverbindungen.

Research for infection protection.





# BODEDEX FORTE

## CHEMICAL COMPOSITION

- Non-ionic and amphoteric surfactants, solvents, complexing agents, corrosion inhibitors, anti-foaming agents, preservatives
- Fragrance, colorant and foam free
- Neutral pH value of diluted solution
- CE Labeling in accordance with the Medical Device Directive

## MATERIAL COMPATIBILITY

- Excellent material compatibility on glass, ceramic, stainless steel, non-ferrous metals, aluminum, hard rubber, hard plastic, silicone, Makrolone, & Plexiglas.

Endorsed by

**FUJINON**  
FUJIFILM

**OLYMPUS**

**STORZ**  
KARL STORZ – ENDOSKOPE

**RICHARD WOLF**

## Bodedex<sup>®</sup> forte

pH-neutral cleaner for heat sensitive and heat resistant instruments. Patented self-acting liquid cleaner

### Product properties

- self acting, strong cleaning effect even with stubborn soiling
- dissolves contrast media residues and biofilms
- special material protectors
- easy to rinse off without leaving residues
- little foaming
- fragrance and colourant free

### Composition

Non-ionic and amphoteric surfactants, solvents, complexing agents, corrosion inhibitors, anti-foaming agents, preservatives

### Directions for use

#### Immersion bath procedure

**Bodedex<sup>®</sup> forte** is used for the pre-cleaning and cleaning of instruments (even of particularly narrow lumen) to optimize the subsequent disinfection and sterilization success. Disassemble or open the instruments respectively and immerge them in **Bodedex<sup>®</sup> forte** solution. The channels have to be filled without any bubbles. Make sure the instruments are completely covered by the solution. After cleaning, rinse all parts with water of at least drinking quality and forward the instruments to further reprocessing.

#### Circulation procedure

**Bodedex<sup>®</sup> forte** can be used in all common semi-automatic circulation procedures for the cleaning process.

#### Ultrasonic bath

**Bodedex<sup>®</sup> forte** can be used in all stainless steel ultrasonic baths.

#### Standing time

The solution should be renewed daily.

### Areas of Application

**Bodedex<sup>®</sup> forte** can be used for self acting cleaning of instruments and equipment in hospitals, primary healthcare and the laboratory.

#### Removal of contrast media residues

Instruments and devices like duodenoscopes that are contaminated with residues of contrast media.

#### Anaesthesia and intensive care

Patient related items such as pharyngeal tubes, masks, laryngoscope spatulas, etc.

#### Surgery

Instruments, surgical equipment / tools, blood aspirators, etc.

#### Endoscopes

Gastrosopes, duodenoscopes, colonoscopes, bronchoscopes, biopsy snares and attachments, etc.

#### Laboratories

Glass and porcelain containers, devices, tubing systems, etc.

#### Urology

Cold light endoscopes (optics) and accessories, rectoscopes, etc.

#### Ophthalmologic instruments

### Listing

CE labelling in accordance with the Medical Device Directive (MDD)

1 % solution :

pH-value approx. 7

### Stability

Diluted solution is stable for one day.

### Characteristics

**Bodedex<sup>®</sup> forte** is a modern, powerful instrument cleaner which is exceptionally suitable for cleaning of medicinal heat-sensitive and heat-resistant instruments as well as laboratory apparatus.

**Bodedex<sup>®</sup> forte** has a neutral pH-value.

Contains also special system of complexing agents, which compensates **water hardness**.

Reliable and gentle removal of stubborn residues of blood, proteins, secretions, grease, biofilms and x-ray contrast media.

**Bodedex<sup>®</sup> forte** assists superior disinfection of instruments and endoscopes.

### Dosage

degree of contamination	dilution of <b>Bodedex<sup>®</sup> forte</b> in water (drinking quality)	contact time
mild	0.5%	5 mins.
heavy	1.0%	10 mins.

Endorsed by



Packing

= 20 bottles / Carton.



Comparison of success when cleaning duodenoscopes (in %)



When cleaning duodenoscopes with a competitive product, 69 % of the devices are free from contamination after complete reprocessing. Thanks to its special formula, Bodedex<sup>®</sup> forte yields a success rate of 98 %.

Research for infection protection.

### Removing biofilm

Organic contamination and biofilm in particular pose a challenge when cleaning heat-sensitive endoscopes. On the one hand, the stubborn soils require high cleaning power; on the other hand, chemicals or too strong mechanical abrasion must not damage the endoscope's sensitive material.

A study shows that, for cleaning flexible endoscopes, the efficacy of a cleaner with neutral pH is superior to enzymatic products (1).

Additionally, according to surveys, the neutral cleaner reduces biofilm with *Pseudomonas aeruginosa* by 1.6 log<sub>10</sub> steps at 40 °C and an exposure time of 5 minutes only, thus being assessed as highly effective for practice (2).

The cleaning success is based on a patented mechanism of action that increases the water solubility of deposits

1 Albrecht H. Infektionspotenzial Kontrastmittelreste: Neuartiges Wirksystem ist bei Duodenoskopen enzymatischen Reinigern überlegen. ENDO-Praxis; 2000; 2:18-20.  
2 Gutachten zur Reinigungswirkung des Neutralreinigers Bodedex forte im manuellen Tauchverfahren gegenüber Biofilm von *Pseudomonas aeruginosa*. HSK, Dr. Horst-Schmidt-Kliniken GmbH, Hygiene-Institut, Wiesbaden, 29.06.2000.



Exclusively Manufactured & Marketed by:  
**RAMAN & WEIL PVT. LTD.**

Registered Office : 15, Chateau Marine, B-Rd., Churchgate, Mumbai - 20.  
Administrative Office : Room No. 36/37, 3rd Floor, Mittal Chamber, Opp. Inox Theatre, Nariman Point, Mumbai - 400 021, India  
Phone : 22049527, 22856397 Email : sales@ramanweil.com  
Visit : www.ramanweil.com



In technical collaboration with:  
**BODE Chemie GmbH**  
Melanchthonstraße 27  
D-22525 Hamburg, Germany  
Must Visit : www.bode-science-center.com  
www.sterillium.com  
Recipient : ISO 9001, ISO 14001, ISO 13485  
EC Certificate of Conformity, WHO - GMP



# BODEDEX FORTE

## COMPARISON OF TWO PRODUCTS

### % of SUCCESS WHEN CLEANING CONTAMINATED DUODENOSCOPES

**69%**

Competitive product (Enzyme based)

**31%**

**98%**

Bodedex Forte

**2**

PUBLISHED STUDY AVAILABLE ON REQUEST

### STANDARD SOLUTION- Use 0.5% to 1%

DEGREE OF CONTAMINATION	DILUTION OF BODEDEX FORTE	CONTACT TIME
MILD	0.5%	5 mins
HEAVY	1.0%	10 mins

### STABILITY

Diluted solution can be used up to 1 day, including high protein loads  
Shelf life of the concentrate is 3 yrs.



500ml

**ALL EFFICACY REPORTS, MATERIAL COMPATIBILITY REPORTS, ENDORSEMENTS AVAILABLE ON REQUEST.**



**Bacillol 25**



500ml, 5L

**AFNOR-NFT -  
72-170**

**Bactericidal  
Fungicidal  
Yeasticidal  
EN 1275  
Viricidal  
COVID-19✓**

**Wettask  
Wipes**



90 Pulls, 30 Pulls

**AFNOR-NFT -  
72-170**

**Bactericidal  
Fungicidal  
Yeasticidal  
EN 1275  
Viricidal  
COVID-19✓**

**Bacttol 25**



225ml, 5L

**Bactericidal  
Fungicidal  
Yeasticidal**

**Gluthyde  
2% &  
2.45%**



5L

**Bactericidal  
Fungicidal  
MRSA  
Sporicidal**

**Gluthyde  
OPA**



5L

**Bactericidal  
Fungicidal  
Tuberculocidal  
Yeasticidal  
Sporicidal**







**HOSPITALS YOU RELY ON,  
RELY ON US**







**AIIMS**

All India Institute Of Medical Sciences



**PGIMER**

Postgraduate Institute of Medical Education and Research, Chandigarh



**CMC Vellore**

Christian Medical College, Vellore



**JIPMER**

Jawaharlal Institute of Post Graduate Medical Education & Research, Puducherry



**COLUMBIA ASIA HOSPITALS**



**SRI RAMACHANDRA**

INSTITUTE OF HIGHER EDUCATION AND RESEARCH

(Category - I Deemed to be University) Porur, Chennai



Lilavati Hospital & Research Centre



**cloudnine**

India's leading chain of maternity hospitals



Catholic Health Association of India

**RELIANCE HOSPITAL** | Every Life Matters



**CARE HOSPITALS**



**Aster**



Breach Candy Hospital Trust



Grant Medical Foundation  
**Ruby Hall Clinic**



**JEHANGIR HOSPITAL**



**STAR HOSPITALS**  
Expertise • Excellence • Empathy



Rainbow Children's Hospital

**Manipal Hospitals**  
LIFE'S ON



**GLOBAL HOSPITALS**  
more to life







Germany

# TRUST IN US

