

MEDICAL GAS SAFETY AND MANAGEMENT



Prepared By :


Mr.BILAL
Mr. MUTHU
Mr. VIJENTHIRAN
Dr. CHANDHANA
Dr.S.AMUDHAGANESH
Ms.LILLY PUSHPAM
Ms.SATHIYA
Dr.VENUGOPAL

Reviewed By :

Dr. Lallu Joseph
Dr. Joseph Fidelis
Dr.Nithya
Dr.Parivalavan Rajavelu

Regulations

- Regulations for use, storage and handling will be according to the authority having jurisdiction, or AHJ
- In the absence of codes, the following may provide guidance:
 - Compressed Gas Association
 - National Fire Protection Association (NFPA)
 - Safety Data Sheet (formerly Material Safety Data Sheet)



Gas Properties

Gases can be:

- Flammable
- Non-flammable
- Oxidizers
- Corrosive
- Asphyxiants
- Poison
- Inert
- Or a mixture



Ignition Temperature

- Ignition temperature: Unique to various solids, vapors and gases, the requisite heat from an open flame source required to ignite materials
- Autoignition temperature: The temperature required to ignite materials absent an open flame source

Gas Containers

- Lecture Bottles
- Cylinders
- Tank Trucks



Gas Containers

- Railroad Tank Cars
- Portable Tanks
- Fixed Storage
- Pipelines




Cylinders

- Construction Must be compatible with the material contained



- Markings
Labeling required to identify the gas in storage and during shipment




Storage Pressure

<u>Types</u>	<u>Storage Pressure (PSI)</u>	<u>Ignition Temperature</u>
•Methane	up to 6000psi	999 °F
•Ethane	544	959 °F
•Propane	109.7	871 °F
•Butane	31	761 °F
•Nitrogen	2,000/below 200 as cryogen	Inert
•Oxygen	2,000/below 200 as cryogen	Inert
•Arsine	219.7	(*see note)




















*Note: Arsine has no given Ignition Temperature but decomposes into arsenic and hydrogen between 446 °F to 464 °F

Color Codes

- Cylinder shells can also be color coded to better identify the contents permitted into the specific type of cylinder
- This eliminates cross-contamination by introducing non-compatible gases into non-specification cylinders



Medical Gas Color Codes

Gas	U.S. Color Code	ISO Color Code
Carbon Dioxide	 Grey	 Grey
He-O₂	 Brown & Green	 Brown & White
Instrument Air	 Red (USA Only)	
Medical Air	 Yellow	 Black & White
Nitrogen	 Black	 Black
Nitrous Oxide	 Blue	 Blue
O₂-He	 Green & Brown	 White & Brown
Oxygen	 Green	 White
Vacuum (Suction)	 White	 Yellow
WAGD (Evac)	 Purple	 Purple

Labels

FTSC Code

Standard numerical code for a gas indicating:


- Flammability
- Toxicity
- State of the gas
- Corrosiveness

CGA V-7 pamphlet provides more in-depth information



Pressure Relief Valve (PRV)

- May be activated by pressure, temperature or spring to permit container contents to escape, thereby averting a container rupture
- The PRV is in the product line



Cylinder Hazards

- Material Hazards

- Flammability
- Spontaneously flammable (arsine, silane and phosphine)
- Corrosivity
- Reactivity
- Poison
- Carcinogenic

- Container Behavior

- Frostbite
- Rupture
- Rocketing
- Boiling liquid expanding vapor explosion (BLEVE)



Oxygen

- Not flammable
- Sensitizes flammable and combustible materials requiring less input heat for ignition
- In some cases, materials impregnated with oxygen can be ignited with static electricity



TYPES OF TERMINAL GAS OUTLETS

Outlets can be installed as flush fitting units, surface-fitting units on booms or pendants, suspended on a hose and gang mounted.

Medical Hanger



Hanging From Ceiling Type Outlet



Ceiling column & pendent



Bedhead Unit



Embedded Type Outlet



Wall Mounted Type Outlet




Plug




Trolley Stand







Medical Cylinder Sizes



GWSCO.COM
(516)334-8200

Cylinder size

	C	D	E	F	G	J
Dimensions (mm)	356×89	457×102	788×102	865×140	1248×40	1450×29
Capacities (L)						
Oxygen	170	340	680	1360	3400	6800
Nitrous oxide	450	900	1800	3600	9000	-
Entonox	-	500	-	2000	5000	-
Air	-	-	-	-	3200	6400
Carbon dioxide	450	1800	1800	-	-	-

Safe Handling & Storage

- Determine safe handling and storage needs based on your industry and the gases with which you work
- Create or follow check lists to best ensure a continuous safety program



Proper Handling

- Use proper hand trucks - do not roll the cylinder on its side
- Provide a forklift cylinder change-out area which maximizes safety for the operator and other staff
- Provide:
 - Ventilation
 - Fire extinguisher
 - PPE



Handling

- Take time to plan what you are going to do with a cylinder and how you are going to do it
- **Always** decide on the side of personal safety



This realistic accident scenario explains why safety shoes shall be worn at cylinder handling.



This accident situation demonstrates, that it is wise to use gloves at cylinder handling

Storage

- Proper ventilation
- Out of the weather
- Not subject to temperature extremes
- Segregate gas types to eliminate fire or chemical reaction hazards
- Use good house keeping practices
- Post signage



Emergency Response Methods

- An extraction hood used for daily operations may be used to vent escaping gas from a cylinder up through a filter
- Hoods and vents may also be equipped with a “scrubber” to neutralize various gases
- Some poison gases may be “scrubbed” this way



Response

- Do you have a trained team?
- Or will you call specialty responders?
- Will special response equipment be needed?
- Special precautions are required for spontaneously combustible gases such as silane



Response

- Determine if you will handle an event alone or with off-site help
- Pre-plan potential zones of harm should your facility have a release
- Practice safety and be safe in handling, use, storage and response to gas incidents



THANK YOU